

# TECHNICAL SPECIFICATION

RESPIRATORY PROTECTION



## Portacount/AccuFIT Quantitative/Qualitative Face Fit Testing Kit & Spares

Qualitative Face Fit Testing Kit / BPT050-000-000



### DESCRIPTION

FaceFit testing kit includes: 2 Nebulizers (1x Sensitivity, 1x Fit Test), 2 Bottles of solution (1x Sensitivity, 1x Fit Test) - Approx 40-60 tests depending on Persons sensitivity, 1x Hood, 1x Collar, 1x Instruction manual and 10x Test report forms.

1. *Full Kit* - BPT050-000-000
2. *Replacement Nebulizers c/w Spare Caps and Cleaning Pin (1x Sensitivity & 1x Fit Test)* - BPT070-000-000
3. *Replacement Solution Bottles (1x Sensitivity, 1x Fit Test)* - BPT080-000-000  
(Approx 40-60 tests depending on Persons sensitivity)
4. *Replacement Hood* - BPT060-000-000
5. *Face Measuring Guide* - BHU000-000-100
6. *Force® Portacount/AccuFIT Quantitative Fit Test Adaptor* - BPT092-000-000

### OPTIONS

#### Qualitative Face Fit Testing Kit (BPT050-000-000)

- Replacement Nebulisers (BPT070-000-000)
- Replacement Test Solutions (BPT080-000-000)
- Replacement Hood (BPT060-000-000)
- Size Guide (BHU000-000-100)
- Force® Portacount / AccuFIT Quantitative Fit Test Adaptor (BPT092-000-000)

### SAFETY STANDARDS

INDG 479 • ED 6273 • ISO 16975-3 • UNI 11719



JSP Limited, Worsham Mill, Minster Lovell, Oxford OX29 0TA, England.

Tel: +44 (0)1993 826050 / Fax: +44 (0)1993 824411 / Email: sales@jsp-safety.com / Web: www.jsp-safety.com

Do Not Publish without written permission from JSP Limited. Contact copyright@jsp-safety.com. Information correct in the United Kingdom at date of issue and may be subject to change. Any third-party machine and automated translations processed by users may have inaccuracies and should not be relied upon.

Page 1 of 2 / Issued 26 Jan 2022 / Copyright 2018 - 2022. JSP Limited. All Rights Reserved.



## SPECIFICATIONS

Material  
Test Hood - Polypropylene with Clear PVC Panel

Material  
Nebuliser Inserts - Polycarbonate

Material  
Fit Test Solution 55ml - HDPE Container, 0.169%  
Denatonium Benzoate Water Solution.

Material  
Carry Box - Correx

Storage Conditions  
Store Solutions Between +5°C & + 40°C



Product Weight  
full Kit 650g

Material  
Nebuliser - Polycarbonate With Polythene & Man  
made Rubber Fittings

Material  
Sensitivity Solution 55ml - HDPE Container,  
0.0135% Denatonium Benzoate Water Solution.

Material  
Nebulisers & Solution Holder - Foam

Latex free  
No

Storage Conditions  
Relative Humidity Not Greater Than 75%

## DOWNLOADS

 [Image Pack \(ZIP\)](#)

 [Declaration of Conformity \(PDF\)](#)



JSP Limited, Worsham Mill, Minster Lovell, Oxford OX29 0TA, England.

Tel: +44 (0)1993 826050 / Fax: +44 (0)1993 824411 / Email: sales@jgpsafety.com / Web: www.jgpsafety.com

Do Not Publish without written permission from JSP Limited. Contact copyright@jgpsafety.com. Information correct in the United Kingdom at date of issue and may be subject to change. Any third-party machine and automated translations processed by users may have inaccuracies and should not be relied upon.

Page 2 of 2 / Issued 26 Jan 2022 / Copyright 2018 - 2022. JSP Limited. All Rights Reserved.

